

Inpatient Medications Pharmacy Interface Automation (PIA)

Installation Guide



Department of Veterans Affairs

September 2016

Version 1.0

Revision History

Date	Version	Description	Author
4/25/16	0.5	Updated Back out and Rollback sections.	Rachel Walters
1/5/16	0.4	Technical edit.	Rachel Walters
12/31/15	0.3	Added Inbound Configuration.	James Harris
11/24/15	0.2	Added Disaster Recovery section.	Rachel Walters
9/15/15	0.1	Initial draft.	Rachel Walters Tony Nixon Mohamed Anwer

Table of Contents

1. Introduction	4
2. System Requirements	4
2.1. Pre-Installation	4
2.2. Patch Installation	5
2.3. Download and Extract Procedure.....	8
2.4. Database Creation	8
2.5. Installation Scripts.....	8
2.6. Cron Scripts	9
2.7. System Configuration.....	9
2.8. Database Tuning	9
3. Implementation	10
4. Back-out Procedure	10
4.1. Back-out Strategy	10
4.2. Back-out Considerations	10
4.2.1. Load Testing	10
4.2.2. User Acceptance Testing	10
4.3. Back-out Criteria	10
4.4. Back-out Risks	10
4.5. Authority for Back-out.....	10
4.6. Back-out Procedure.....	11
5. Rollback Procedure	16
5.1. Rollback Considerations.....	16
5.2. Rollback Criteria	16
5.3. Rollback Risks	16
5.4. Authority for Rollback	16
5.5. Rollback Procedure	16
6. Disaster Recovery and Continuity of Operations	18

1. Introduction

This document provides installation instructions for Clinical Ancillary Services (CAS) Development Delivery of Pharmacy enhancements (DDPE) Pharmacy Interface Automation (PIA).

2. System Requirements

The following are software elements for Inpatient Medications.

Inpatient Medications Software Elements

Application	Version
Adverse Reaction Tracking	4.0
Decision Support System	3.0
Inpatient Medications	5.0
Kernel	8.0
Laboratory	5.2
Mailman	8.0
National Drug File	4.0
Nursing	4.0
Order Entry/Results Reporting	3.0
Outpatient Pharmacy	7.0
Patient Information Management Systems	5.3
Pharmacy Data Management	1.0
RPC Broker (32-bit)	1.1
Toolkit	7.3
VA FileMan	22.0

2.1. Pre-Installation

The following two patches are provided for this project and should be installed in the order listed below, as the logical link, PSJ PADE that is used for this project is transported in the PSJ patch.

PSJ*5*317 - Inpatient Medications V. 5.0

PSS*1*193 - Pharmacy Data Management V. 1.0

The associated patches for PSJ*5*317 are the following:

PSJ*5*191

PSJ*5*228

PSJ*5*244

PSJ*5*281

PSJ*5*285

The associated patch for PSS*1*193 is PSS*1*180.

2.2. Patch Installation

The scope of the Pharmacy Interface Automation Project (PIA) will provide a new standard bi-directional HL7 interface to the Pharmacy Automated Dispensing Equipment (PADE) located at the point of care areas such as Inpatient wards, Outpatient Clinics etc. Depending on the PADE setup, outbound HL7 messages will be sent to the respective PADE via the VistA Interface Engine (VIE), when the following event happens:

- In Inpatient Medications, when orders are verified, edited, renewed, discontinued, reinstated, put on hold/removed from hold etc.
- In Computerized Patient Record System (CPRS), when orders are discontinued.
- In the REGISTRATION package when patients are admitted, discharged, transferred, cancel admit, cancel discharge, cancel transfer, bed switch etc.
- In the Scheduling package when patients are checked in.
- In Pharmacy Data Management, when new drugs are entered or modified.
- Options to schedule nightly jobs to run daily to send patient clinic appointments or surgery cases.
- Options to send medication orders for a single patient or by clinics or by wards.

Installation Instructions

PSJ*5*317 - Inpatient Medications V. 5.0.

This patch should be installed with users off the system during off-peak hours. Installation takes less than five minutes.

1. Use the INSTALL/CHECK MESSAGE option on the PackMan menu.
2. From the Kernel Installation & Distribution System menu, select the Installation menu.
3. From this menu, you may select to use the following options:
(when prompted for INSTALL NAME, enter PSJ*5.0*317)
 - a. Backup a Transport Global - this option will create a backup message of any routines exported with the patch. It will NOT backup any other changes such as DDs or templates.
 - b. Compare Transport Global to Current System - this option will allow you to view all changes that will be made when the patch is installed. It compares all components of the patch (routines, DDs, templates, etc.).
 - c. Verify Checksums in Transport Global - this option will ensure

the integrity of the routines that are in the transport global.

4. Use the Install Package(s) option and select the package PSJ*5.0*317.
5. When prompted "Want KIDS to Rebuild Menu Trees Upon Completion of Install? NO//" respond NO.
6. When prompted "Want KIDS to INHIBIT LOGONs during the install? NO//" respond NO.
7. When prompted "Want to DISABLE Scheduled Options, Menu Options, and Protocols? NO//" respond NO.

Routine Information:

=====

The second line of each of these routines now looks like:

```
;5.0;INPATIENT MEDICATIONS; **[Patch List]**;16 DEC 97;Build 125
```

The checksums below are new checksums, and
can be checked with CHECK1^XTSUMBLD.

Routine Name: PSGBRJ		
Before: B19761653	After: B21011842	**12,50,244,317**
Routine Name: PSGOE7		
Before: B30726260	After: B40214304	**9,26,34,52,55,50,87,111,181, 254,267,260,288,281,317**
Routine Name: PSGOE82		
Before: B21878874	After: B33061347	**2,35,50,67,58,81,127,168,181, 276,317**
Routine Name: PSGOE92		
Before: B30935948	After: B43137679	**2,35,50,58,81,110,215,237, 276,316,317**
Routine Name: PSGOEF1		
Before: B28843982	After: B36069326	**2,7,35,39,45,47,50,63,67,58, 95,110,186,181,267,317**
Routine Name: PSGPEN		
Before: B39805137	After: B57109712	**30,37,50,58,115,110,127,129,317**
Routine Name: PSGPLR		
Before: B39066813	After: B40061443	**10,50,67,119,129,191,317**
Routine Name: PSJ317P		
Before: n/a	After: B6017554	**317**
Routine Name: PSJHLU		
Before: B45959307	After: B50964072	**1,56,72,102,134,181,267,285,317**
Routine Name: PSJLMPRU		
Before: B14898207	After: B19398098	**16,58,85,110,185,181,267,317**
Routine Name: PSJLMUDE		
Before: B66391198	After: B84078277	**7,47,50,63,64,58,80,116,110, 111,164,175,201,181,254,267, 228,317**
Routine Name: PSJO		

Before: B28881243	After: B32312057	**31,58,110,181,267,275,317**
Routine Name: PSJ02		
Before: B20134007	After: B21667076	**58,317**
Routine Name: PSJPAD50		
Before: n/a	After: B143814904	**317**
Routine Name: PSJPAD70		
Before: n/a	After: B193106268	**317**
Routine Name: PSJPAD7I		
Before: n/a	After: B91989090	**317**
Routine Name: PSJPAD7U		
Before: n/a	After: B183953635	**317**
Routine Name: PSJPADE		
Before: n/a	After: B87815402	**317**
Routine Name: PSJPADIT		
Before: n/a	After: B213276051	**317**
Routine Name: PSJPADPT		
Before: n/a	After: B74411481	**317**
Routine Name: PSJPADSI		
Before: n/a	After: B207002155	**317**
Routine Name: PSJPDAPP		
Before: n/a	After: B26391444	**317**
Routine Name: PSJPDCL		
Before: n/a	After: B58778000	**317**
Routine Name: PSJPDCLA		
Before: n/a	After: B123299808	**317**
Routine Name: PSJPDCLU		
Before: n/a	After: B182392496	**317**
Routine Name: PSJPDRIN		
Before: n/a	After: B220627483	**317**
Routine Name: PSJPDRIP		
Before: n/a	After: B92189032	**317**
Routine Name: PSJPDRTTP		
Before: n/a	After: B164565338	**317**
Routine Name: PSJPDRTTR		
Before: n/a	After: B204379013	**317**
Routine Name: PSJPDRTU1		
Before: n/a	After: B195220461	**317**
Routine Name: PSJPDRTUT		
Before: n/a	After: B233191952	**317**

PSS*1*193 – Pharmacy Data Management V. 1.0.

This patch should be installed with users off the system during off-peak hours. Installation takes less than two minutes.

1. Use the INSTALL/CHECK MESSAGE option on the PackMan menu.
2. From the Kernel Installation & Distribution System menu, select the Installation menu.
3. From this menu, you may select to use the following options:
(when prompted for INSTALL NAME, enter PSS*1.0*193)

- a. Backup a Transport Global - this option will create a backup message of any routines exported with the patch. It will NOT backup any other changes such as DDs or templates.
 - b. Compare Transport Global to Current System - this option will allow you to view all changes that will be made when the patch is installed. It compares all components of the patch (routines, DDs, templates, etc.).
 - c. Verify Checksums in Transport Global - this option will ensure the integrity of the routines that are in the transport global.
4. Use the Install Package(s) option and select the package PSS*1.0*193.
 5. When prompted "Want KIDS to INHIBIT LOGONs during the install? NO//" respond NO.
 6. When prompted "Want to DISABLE Scheduled Options, Menu Options, and Protocols? NO//" respond NO.

Routine Information: =====

The second line of each of these routines now looks like:

```
;;1.0;PHARMACY DATA MANAGEMENT; **[Patch List]**;9/30/97;Build 17
```

The checksums below are new checksums, and can be checked with CHECK1^XTSUMBLD.

Routine Name: PSSDEE

```
Before: B98607664 After: B110481507 **3,5,15,16,20,22,28,32,34,33,
38,57,47,68,61,82,90,110,155,
156,180,193**
```

Routine Name: PSSHLDFS

```
Before: n/a After: B33701533 **193**
```

Routine Name: PSSMSTR

```
Before: B1853023 After: B51317382 **82,193**
```

2.3. Download and Extract Procedure

N/A

2.4. Database Creation

N/A

2.5. Installation Scripts

N/A

2.6. Cron Scripts

N/A

2.7. System Configuration

2.7.1. Outbound

Outbound to PADE setup will be configured in VistA to send the information to the appropriate PADE. The main purpose of the Outbound setup is to map the send areas to the location of the cabinets in the wards, clinics, and/or operating rooms.

Please refer to the [Pharmacy Interface Automation Startup and Troubleshooting Guide](#) for more details.

2.7.2. Inbound

The PADE Inbound HL7 interface must be configured prior to first use. The Inbound HL7 interface does not require the PADE Outbound HL7 interface. However, the inbound interface does require an entry in the outbound system file PADE SYSTEM SETUP (#58.7) before it can be activated.

Please refer to the [Pharmacy Interface Automation Startup and Troubleshooting Guide](#) for more details.

2.8. Database Tuning

N/A

3. Implementation

The PIA Implementation Plan specifies how the Pharmacy Interface Automation project will be evaluated and subsequently deployed within the Office of Enterprise Development (OED).

The PIA Implementation Plan:

- a. Describes the phased implementation for the implementation
- b. States objectives for the implementation
- c. Addresses risks associated with implementation

Refer to the [PIA Implementation Plan](#) for more detailed information.

4. Back-out Procedure

The following section details the back-out procedure for PIA.

4.1. Back-out Strategy

See section 4.6 for more details.

4.2. Back-out Considerations

Back-out pertains to a return to the last known good operational state of the software and appropriate platform settings.

4.2.1. Load Testing

To be determined

4.2.2. User Acceptance Testing

User Acceptance Testing is in progress.

4.3. Back-out Criteria

The project is canceled and the implemented features are no longer wanted by the stake holders.

4.4. Back-out Risks

To be determined

4.5. Authority for Back-out

The authority would come from the IPT and the VA project manager.

4.6. Back-out Procedure

In the event that the Pharmacy Interface Automation enhancements must be backed out, the modified routines must have been backed-up during patch installation using the following option:

Backup a Transport Global [XPD BACKUP]

This option creates a MailMan message of any routines exported with this patch. (If you need to preserve components that are not routines, you must back them up separately.)

Select Kernel Installation & Distribution System Option: Installation

- 1 Load a Distribution
- 2 Verify Checksums in Transport Global
- 3 Print Transport Global
- 4 Compare Transport Global to Current System
- 5 Backup a Transport Global**
- 6 Install Package(s)
- Restart Install of Package(s)
- Unload a Distribution

Select Installation <TEST ACCOUNT> Option: **5 Backup a Transport Global**

select INSTALL NAME: **Backup of Patch_XXXXX** 12/5/15@13:29:01

=> **Backup of Patch_XXXXX**

This Distribution was loaded on Feb 05, 2013@13:29:01 with header of
Patch_XXXXX Test Version

It consisted of the following Install(s): **Patch_XXXXX Test Version**

Subject: **Backup of Patch_XXXXX** install on Feb 05, 2013

Replace

Loading Routines for XXXXX

Routine YYYY1

Routine YYYY2

Routine YYYY3

Send mail to: ADPAC,ONE// ADPAC,ONE

Select basket to send to: IN// BACKUP PATCH (Folder for FORUM)

And Send to:

Restore Pre-Patch Routines (MailMan)

Go to the **Backup of Patch_XXXXX** message in Mailman.

At the Enter message action prompt, enter “X” to “Xtract PackMan”

At the Select PackMan Function prompt enter the number **6** to **Install/Check Message**

At the end of this process the pre-patch routines are restored.

Note: See header “Install the Patch Backup” for detail

Install the Patch Backup

BACKUP PATCH Basket, 144 messages (1-144), 117 new

*=New/!=Priority.....Subject.....From...
141. **Backup of Patch_XXXXX** install on Oct 01, 2015 **ADPAC,ONE**

IN Basket Message: 1// 141

Subj: **Backup of XXXXX** install on Oct 01, 2015 [#000000] 10/01/15@11:14
2016 lines

From: **ADPAC,ONE** In 'IN' basket. Page 1

\$TXT PACKMAN BACKUP Created on Thursday, 10/1/15 at 11:14:50 by ADPAC,ONE

Enter message action (in IN basket): Ignore// **Xtract PackMan**
Select PackMan function: **6 INSTALL/CHECK MESSAGE**

Warning: Installing this message will cause a permanent update of globals
and routines.

Do you really want to do this? NO// **YES**

Routines are the only parts that are backed up. NO other parts
are backed up, not even globals. You may use the 'Summarize Message'
option of PackMan to see what parts the message contains.
Those parts that are not routines should be backed up separately
if they need to be preserved.

Shall I preserve the routines on disk in a separate back-up message? YES// **NO**
No backup message built.

Line 123 Message #000000 Unloading Routine Routine_Name1 (PACKMAN_BACKUP)
Line 345 Message #000000 Unloading Routine Routine_Name2 (PACKMAN_BACKUP)
Line 567 Message #000000 Unloading Routine Routine_Name3 (PACKMAN_BACKUP)
Line 789 Message #000000 Unloading Routine Routine_Name4 (PACKMAN_BACKUP)
=====

New Routine(s)

New routines implemented by the patches can be deleted/removed by using the following option:
Delete Routines [XTRDEL]

This option can be found under the **Routine Tools** menu

Select OPTION NAME: **XUPROG** **Programmer Options**

KIDS Kernel Installation & Distribution System ...

```

NTEG   Build an 'NTEG' routine for a package
PG      Programmer mode
        Calculate and Show Checksum Values
        Delete Unreferenced Options
        Error Processing ...
        Global Block Count
        List Global
        Map Pointer Relations
        Number base changer
        Routine Tools ...
        Test an option not in your menu
        Verifier Tools Menu ...

```

Select Programmer Options <TEST ACCOUNT> Option: **ROUTINE Tools**

```

%Index of Routines
Check Routines on Other CPUs
Compare local/national checksums report
Compare routines on tape to disk
Compare two routines
Delete Routines
First Line Routine Print
Flow Chart Entire Routine
Flow Chart from Entry Point
Group Routine Edit
Input routines
List Routines
Load/refresh checksum values into ROUTINE file
Output routines
Routine Edit
Routines by Patch Number
Variable changer
Version Number Update

```

Select Routine Tools <TEST ACCOUNT> Option: **DELETE Routines**

```

ROUTINE DELETE
All Routines? No => No
Routine: ROUT999
Routine:
1 routine
1 routines to DELETE, OK: NO// YES
ROUT999
Done.

```

Other Components

Data dictionary and template modifications must be removed using a follow-up patch.

<u>File Name (#)</u>	<u>Field Name (#)</u>	<u>New/Modified/Deleted</u>
PADE SYSTEM SETUP (#58.7)	All fields	New
PADE SEND AREA (#58.71)	All fields	New
PADE OUTBOUND MESSAGES (#58.72)	All fields	New
PADE DISPENSING DEVICE (#58.63)	All fields	New

PADE INBOUND TRANSACTIONS (#58.6)	All fields	New
PADE INVENTORY SYSTEM (#58.601)	All fields	New
PADE USER (#58.64)	All fields	New
INPATIENT WARD PARAMETERS (#59.6)	DEFAULT 0 ON PADE	New
	PRE-EXCHANGE (#8)	

Kernel Parameters	<u>New/Modified/Deleted</u>
PSJ PADE OE BALANCES	New

Mail Group	<u>New/Modified/Deleted</u>
PSJ PADE DISPENSE ALERTS	New

Options	<u>New/Modified/Deleted</u>
Unit Dose Medications [PSJU MGR]	Modified Menu
PADE Main Menu [PSJ PADE MAIN MENU]	New Menu
PADE Send Area Setup [PSJ PADE SEND AREA SETUP]	New Option
PADE System Setup [PSJ PADE SETUP]	New Option
PADE Inventory Setup [PSJ PADE INVENTORY MENU]	New Menu
Inventory System Setup [PSJ PADE INVENTORY SYSTEM]	New Option
Dispensing Device Setup [PSJ PADE DEVICE SETUP]	New Option
PADE Send Surgery Cases [PSJ PADE SEND SURGERY CASES]	New Option
PADE Surgery Task [PSJ PADE SURGERY TASK]	New Option
PADE Reports [PSJ PADE REPORTS MENU]	New Menu
PADE On-Hand Amounts [PSJ PADE INVENTORY REPORT]	New Option
PADE Transaction Report [PSJ PADE TRANSACTION REPORT]	New Option
PADE System Division Setup [PSJ PADE DIVISION SETUP]	New Option
PADE Send Patient Orders [PSJ PADE SEND ORDERS]	New Option
PSJ PADE Appointment Task [PSJ PADE APPOINTMENT TASK]	New Option
Send Drug File Entries to External Interface [PSS MASTER FILE ALL]	Modified Option

Protocols	<u>New/Modified/Deleted</u>
PSJ ADT-A01 CLIENT	New
PSJ ADT-A01 ROUTER	New
PSJ ADT-A01 SERVER	New
PSJ ADT-A02 SERVER	New
PSJ ADT-A02 CLIENT	New
PSJ ADT-A03 SERVER	New
PSJ ADT-A03 CLIENT	New
PSJ ADT-A11 SERVER	New
PSJ ADT-A11 CLIENT	New
PSJ ADT-A12 SERVER	New
PSJ ADT-A12 CLIENT	New
PSJ ADT-A13 SERVER	New
PSJ ADT-A13 CLIENT	New
PSJ SIU-SDAM ROUTER	New
PSJ SIU-S12 SERVER	New
PSJ SIU-S12 CLIENT	New
PSJ RDE011 SERVER	New
PSJ RDE011 CLIENT	New
PSJ PADE OMS-005 EVENT	New

PSJ PADE OMS-005 SUB	New
PSJ PADE OMS-005 EVENT 2.3	New
PSJ PADE OMS-005 SUB 2.3	New
PSJ LM PADE ACTIVITY	New
PSJ LM PROFILE HIDDEN ACTIONS	Modified
PSS MFNM01 SERVER	New
PSS MFNM01 CLIENT	New

<u>Security Keys</u>	<u>New/Modified/Deleted</u>
-----------------------------	------------------------------------

PSJ PADE ADV	New
PSJ PADE MGR	New
PSS PADE INIT	New

<u>Templates</u>	<u>New/Modified/Deleted</u>
-------------------------	------------------------------------

PSJ PADE SYSTEM	Input	58.7	New
PSJ PADE INVENTORY	Input	58.601	New
PSJ PADE DISPENSING DEVICE	Input	58.63	New

5. Rollback Procedure

The following sections detail the rollback procedure for PIA.

5.1. Rollback Considerations

A follow-up patch would be needed to remove new data entries established by the data dictionary. Specific rollback details will be incorporated in subsequent versions, but a patch would only be created if necessary.

Note: These new data entries are the result of a new field added to the data dictionary or the modification of an existing field in the data dictionary.

5.2. Rollback Criteria

The back-out of PHARMACY INTERFACE AUTOMATION (PIA) patches that modified existing fields, and established new fields would be justification for rollback.

5.3. Rollback Risks

None

5.4. Authority for Rollback

The authority would come from the IPT and the VA project manager.

5.5. Rollback Procedure

A follow-up patch for each namespace would be needed to delete/modify that namespace's data dictionary entries that were added/modified and other non-routine components added/modified by this projects patches and would follow the basic logic flow below.

Logic flow using fileman API calls as much as possible for the below actions

1. New data fields would need to be erased that were likely populated by the new functionality, while the new data fields are still valid in the data dictionary.
2. New file cross references would now need to be deleted.
3. New fields need to be deleted out of the data dictionary.
4. Modified data dictionary fields would need to be restored.
5. New SECURITY KEY file (#19.1) entries needs to be deleted.
6. Existing PROTOCOL file (#101) entries would need to be restored while new file entries would be deleted from the site.
7. Existing OPTION file (#19) entries would need to be restored while new file entries would be deleted from the site.
8. Existing PARAMETERS file (#8989.5) entries would need to be restored while new file entries would be erased and then deleted from the Kernel Parameter Definition file at the site.

9. Existing PARAMETER DEFINITION file (#8989.51) entries would need to be restored while new file entries would be erased and then deleted from the PARAMETER DEFINITION file at the site.
10. New INPUT TEMPLATES file (.402) entries need to be deleted.
11. New HL LOGICAL LINK file (#870) entry needs to be deleted.
12. New HL7 APPLICATION PARAMETER file (#771) entries needs to be deleted.

6. Disaster Recovery and Continuity of Operations

Each VistA facility as well as regional data centers are responsible for their own disaster recovery (DR) and Continuity of Operations (COOP). Please refer to the VistA Disaster Recovery and Continuity of Operations Plans at the specific facility. Most of these documents are considered confidential in that information that could disrupt any of these facilities DR or COOP could cause catastrophic data loss. Therefore the following link is just an example:

- http://vaww.oed.portal.va.gov/projects/ecms_to_ifcap/Phase%201%20Archive/References/VistA%20DRP%20and%20CP/VANIHCS%20VISTA%20DRP%20November%202011.docx
- <http://vaww.eie.va.gov/techstrategy/TAR/SEDR120578/VA%20Enterprise%20Disaster%20Recovery%20Service%20Tiers%20Standard%2008-16-2012%20Final.pdf>