



**Medication Order Check Healthcare Application  
(MOCHA) Enhancements 2  
Combined Build 1.0**

**Installation Guide**

**(Rollback Plan)**

GMRA\*4\*46, OR\*3\*269, PSJ\*5\*281 and PSO\*7\*411  
(MOCHA Enhancements 2 Combined Build 1.0)

PSS\*1\*175 (Stand Alone)

April 2016

---

Department of Veterans Affairs

Product Development

*(This page included for two-sided copying.)*

## Table of Contents

<b>1. Introduction .....</b>	<b>1</b>
1.1 Purpose.....	1
<b>2. Pre-Requisite Considerations.....</b>	<b>1</b>
2.1 Minimum Required Packages .....	1
2.2 Required Patches.....	2
<b>3. Installation Considerations/Restrictions.....</b>	<b>3</b>
<b>4. Installation Order.....</b>	<b>5</b>
<b>5. Rollback Preparation.....</b>	<b>7</b>
<b>6. Installation of PSS*1*175 .....</b>	<b>9</b>
6.1 Installation Steps .....	9
6.1.1 Install Example: .....	9
<b>7. Installation of MOCHA_ENH_2_GMRA_OR_PSI_PSO.KID .....</b>	<b>11</b>
7.1 Pre-Installation Instructions .....	11
7.2 Installation Steps.....	11
7.2.1 Install Example: .....	12
7.3 Routine Deletions.....	16
<b>8. Rollback Procedure.....</b>	<b>17</b>
8.1 Rollback Requirements.....	17
8.1.1 Software Required for Rollback .....	17
8.1.2 Overview of Rollback Build .....	17
8.2 Rollback Considerations/Restrictions.....	18
8.3 Rollback Steps .....	18
<b>Appendix: Acronyms .....</b>	<b>21</b>



# 1. Introduction

Medication Order Check Healthcare Application (MOCHA) Enhancements (Enh) 2 provides additional functionality for the existing MOCHA 1.0 drug interaction and duplicate therapy enhanced order check functionality. In addition, modifications have been made for Clinical Reminder Order Checks (CROC) and allergy functionality.

The following patches make up the release of MOCHA Enh 2:

Patch	Type	Description
PSS_1_175.KID	Host File	MOCHA UFT ENHANCEMENTS 2 PSS - Pharmacy Data Management patch PSS*1*175
MOCHA_ENH_2_GMRA_OR_PSJ_PSO.KID	Host File	MOCHA UFT ENHANCEMENTS 2 combined build that contains the following patches: <ul style="list-style-type: none"><li>• GMRA*4*46 - Adverse Reaction Tracking patch</li><li>• OR*3*269 - Computerized Patient Record System(CPRS) patch</li><li>• PSJ*5*281 - Inpatient Medications patch</li><li>• PSO*7*411 - Outpatient Pharmacy patch</li></ul>

## 1.1 Purpose

The purpose of this Installation Guide is to provide installation steps for MOCHA Enhancements 2. The intended audience for this document is the Information Resources Management Service (IRMS) staff responsible for installing these patches.

## 2. Pre-Requisite Considerations

No pre-requisite considerations.

### 2.1 Minimum Required Packages

The patches described in this installation guide can only be run with a standard Massachusetts General Hospital Utility Multi-Programming System (MUMPS) operating system and require the following Department of Veterans Affairs (VA) software packages.

Package	Minimum Version Needed
Pharmacy Data Management	1.0
VA FileMan	22.0
Kernel	8.0
HealtheVet Web Services Client (HWSC)	1.0
Outpatient Pharmacy	7.0
Inpatient Medications	5.0
Order Entry/Results Reporting	3.0
Adverse Reaction Tracking	4.0

The above software must be installed for these patches to be completely functional.

## 2.2 Required Patches

The following patches should already be installed on your system:

- GMRA\*4\*10
- GMRA\*4\*44
- OR\*3\*286
- OR\*3\*306
- OR\*3\*311
- OR\*3\*346
- OR\*3\*352
- OR\*3\*395
- PSJ\*5\*160
- PSJ\*5\*199
- PSJ\*5\*278
- PSJ\*5\*286
- PSJ\*5\*299
- PSJ\*5\*305
- PSJ\*5\*309
- PSJ\*5\*320
- PSO\*7\*294
- PSO\*7\*313
- PSO\*7\*372
- PSO\*7\*390
- PSO\*7\*416
- PSO\*7\*429
- PSO\*7\*436
- PSO\*7\*438
- PSS\*1\*169

Important: Recently released patch GMRA\*4\*48 contains an Allergy Assessment Clean Up utility tool that produces a report that identifies discrepancies between the ADVERSE REACTION ASSESSMENT file (#120.86) and the PATIENT ALLERGIES file (#120.8). It is imperative that GMRA\*4\*48 be installed prior to the installation of these MOCHA Enhancements 2 patches, and that all data discrepancies identified by GMRA\*4\*48 have been corrected. If not, then installing the MOCHA Enhancements 2 patches may result in drug allergy order checks not displaying in CPRS and pharmacy during the medication order entry processes.

### 3. Installation Considerations/Restrictions

These patches should be installed when Pharmacy applications are not in use, no other pharmacy patches are being installed and when tasked jobs from Clinical Applications are not running.

Installation should also occur when CPRS usage is at a minimum, particularly medication activities and when allergy data is not being entered or edited in CPRS.

The length of time for the install is dependent upon the number of entries in the ORDER CHECK INSTANCES (#100.05) file that will be converted and could range from five to thirty minutes.

**Do not queue the installation of this patch.**

Failure to heed this warning may result in 'source routine edited' errors during a database update. Edits may be lost and records may be left in an inconsistent state.





## 4. Installation Order

MOCHA Enh 2 patches should be installed in the following order and should be installed consecutively together. Please review individual patch install information in the following sections for pre and post installation instructions.

1. PSS\_1\_175.KID – HOST file with patch name:

PSS\*1.0\*175

2. MOCHA\_ENH\_2\_GMRA\_OR\_PSO.KID – HOST file with patch name:

MOCHA Enh 2 COMBINED BUILD 1.0



## 5. Rollback Preparation

The rollback procedure for these patches should only occur when there is concurrence from the Enterprise Product Support and MOCHA development teams, because of the complexity and risk involved in a rollback. Normal installation back-ups using KIDS will back up only Mumps routines. For all non-routine components of these builds, Enterprise Product Support will have a build available if needed. Make sure the 'Backup a Transport Global' step in sections 6 and 7 of this document is followed, so you do have a backup of all the routines if needed.

If your site has any local modifications to any of the following non-routine components of the builds, they will need to be backed up locally. Please contact Enterprise Product Support prior to installing if there are questions.

- ORDER CHECK INSTANCES (#100.05) File Data Dictionary
- PSJ DISPLAY DRUG ALLERGIES Protocol




## 6. Installation of PSS\*1\*175

### 6.1 Installation Steps

1. Download PSS\_1\_175.KID into your local directory.
2. From the Kernel Installation and Distribution System (KIDS) Menu, select the Installation menu.
3. Use **Load a Distribution**. You may need to prepend a directory name. When prompted for "Enter a Host File:", respond with PSS\_1\_175.KID.

Example: USER\$:[ABC]PSS\_1\_175.KID

4. From this menu, you may then select to use the following options:

<b>Note:</b>	
	The following are OPTIONAL - (When prompted for the INSTALL NAME, enter PSS*1.0*175).

- **Verify Checksums in Transport Global** - This option will allow you to ensure the integrity of the routines that are in the transport global.
  - **Print Transport Global** - This option will allow you to view the components of the KIDS build.
  - **Compare Transport Global to Current System** - This option will allow you to view all changes that will be made when this patch is installed. It compares all components of this patch (routines, DD's, templates, etc.).
5. Select the installation option **Backup a Transport Global**. This option will create a backup message of any routines exported with this patch in case you need to backout this patch. It will not backup any other changes such as Data Dictionaries (DD's) or templates. It is important this step be followed, because if backout of this patch is necessary, having this backup will make the process much easier.
  6. Select the Installation option **Install Package(s)**. This is the step to start the installation of this KIDS patch:
    - a. Choose the Install Package(s) option to start the patch install and enter PSS\*1.0\*175 at the INSTALL NAME prompt.
    - b. When prompted 'Want KIDS to INHIBIT LOGONs during the install? NO//' answer NO.
    - c. When prompted 'Want to DISABLE Scheduled Options, Menu Options, and Protocols? NO//' answer NO.

#### 6.1.1 Install Example:

```
Select Installation <TEST ACCOUNT> Option: 6  Install Package(s)
Select INSTALL NAME:      PSS*1.0*175      5/20/14@15:13:40
=> PSS*1*175

This Distribution was loaded on May 20, 2014@15:13:40 with header of
```

```

PSS*1*175
It consisted of the following Install(s):
PSS*1.0*175
Checking Install for Package PSS*1.0*175

Install Questions for PSS*1.0*175
Want KIDS to INHIBIT LOGONs during the install? NO//
Want to DISABLE Scheduled Options, Menu Options, and Protocols? NO//

Enter the Device you want to print the Install messages.
You can queue the install by enter a 'Q' at the device prompt.
Enter a '^' to abort the install.

DEVICE: HOME//    SSH VIRTUAL TERMINAL

                                PSS*1.0*175

Install Started for PSS*1.0*175 :
                        May 20, 2014@15:26:49

Build Distribution Date: May 20, 2014

Installing Routines:
                        May 20, 2014@15:26:49

Updating Routine file...

Updating KIDS files...

PSS*1.0*175 Installed.
                        May 20, 2014@15:26:50

Not a production UCI

NO Install Message sent

Install Completed

```

## 7. Installation of MOCHA\_ENH\_2\_GMRA\_OR\_PSJ\_PSO.KID

### 7.1 Pre-Installation Instructions


PSS\*1.0\*175 must have been successfully installed prior to installation of this host file.

### 7.2 Installation Steps

1. Download MOCHA\_ENH\_2\_GMRA\_OR\_PSJ\_PSO.KID into your local directory.
2. From the Kernel Installation and Distribution System (KIDS) Menu, select the Installation menu.
3. Use **Load a Distribution**. You may need to prepend a directory name. When prompted for "Enter a Host File:", respond with MOCHA\_ENH\_2\_GMRA\_OR\_PSJ\_PSO.KID.

Example: USER\$:[ABC]MOCHA\_ENH\_2\_GMRA\_OR\_PSJ\_PSO.KID

4. From this menu, you may then select to use the following options:

<b>Note:</b>	
	The following are OPTIONAL - (When prompted for the INSTALL NAME, enter MOCHA ENH 2 COMBINED BUILD 1.0).

- **Verify Checksums in Transport Global** - This option will allow you to ensure the integrity of the routines that are in the transport global.
  - **Print Transport Global** - This option will allow you to view the components of the KIDS build.
  - **Compare Transport Global to Current System** - This option will allow you to view all changes that will be made when this patch is installed. It compares all components of this patch (routines, DD's, templates, etc.).
5. Select the installation option Backup a Transport Global. This option will create a backup message of any routines exported with this patch in case you need to backout this patch. It will not backup any other changes such as Data Dictionaries (DD's) or templates. It is important this step be followed, because if backout of this patch is necessary, having this backup will make the process much easier.
  6. Select the Installation option **Install Package(s)**. This is the step to start the installation of this KIDS patch:
    - a. Choose the Install Package(s) option to start the patch install and enter MOCHA ENH 2 COMBINED BUILD 1.0 at the INSTALL NAME prompt.
    - b. When prompted 'Want KIDS to Rebuild Menu Trees Upon Completion of Install? NO//', answer NO.
    - c. When prompted 'Want KIDS to INHIBIT LOGONs during the install? NO//' answer NO.

- d. When prompted 'Want to DISABLE Scheduled Options, Menu Options, and Protocols? NO//' answer YES.
- e. When prompted 'Enter options you wish to mark as 'Out Of Order':' answer GMRA PATIENT A/AR EDIT.
- f. When prompted 'Enter protocols you wish to mark as 'Out Of Order':' press the Enter key.
- g. When prompted 'Delay Install (Minutes): (0-60): 0//' press the Enter key.

Please note that because the installation process of this build will run the CPRS Order Check routine compiler, the patch installer will receive a VistA mail message similar to this mail message indicating the Order Check compiler has run to completion:

```

Subj: Order Check Compiler Status  [#77714] 03/24/16@15:23  14 lines
From: POSTMASTER  In 'IN' basket.  Page 1
-----

The Order Check routine compiler has completed normally
on MAR 24,2016 at 15:23 by [nnnnnn]  INSTALLER,PATCH.

ORDER CHECK EXPERT version 1.01 released OCT 29,1998

                Elapsed time: 0 minutes 7 seconds
                        Automatic Mode
        Execution Trace Mode: OFF
Elapsed time Logging Mode:
        Raw Data Logging Mode: OFF
        Lines of code generated: No longer tracked

```

## 7.2.1 Install Example:

```

Select Installation <TEST ACCOUNT> Option: 6  Install Package(s)
Select INSTALL NAME:      MOCHA ENH 2 COMBINED BUILD 1.0      8/11/15@15:05:30
      => MOCHA ENH 2 COMBINED BUILD (8.11.15)  ;Created on Aug 11, 2015@14:

This Distribution was loaded on Aug 11, 2015@15:05:30 with header of
      MOCHA ENH 2 COMBINED BUILD (8.11.15)  ;Created on Aug 11, 2015@14:14:04
It consisted of the following Install(s):
MOCHA ENH 2 COMBINED BUILD 1.0      GMRA*4.0*46      OR*3.0*269      PSJ*5.0*281
      PSO*7.0*411
Checking Install for Package MOCHA ENH 2 COMBINED BUILD 1.0

Install Questions for MOCHA ENH 2 COMBINED BUILD 1.0

Checking Install for Package GMRA*4.0*46
Install Questions for GMRA*4.0*46

Checking Install for Package OR*3.0*269
Install Questions for OR*3.0*269

Incoming Files:

```



```

100.05    ORDER CHECK INSTANCES
Note:  You already have the 'ORDER CHECK INSTANCES' File.

Checking Install for Package PSJ*5.0*281

Install Questions for PSJ*5.0*281

Checking Install for Package PSO*7.0*411

Install Questions for PSO*7.0*411

Want KIDS to INHIBIT LOGONs during the install? NO//
Want to DISABLE Scheduled Options, Menu Options, and Protocols? NO// YES

Enter options you wish to mark as 'Out Of Order': GMRA PATIENT A/AR EDIT      E
nter/Edit Patient Reaction Data

Enter options you wish to mark as 'Out Of Order':

Enter protocols you wish to mark as 'Out Of Order':

Delay Install (Minutes):  (0-60): 0//

Enter the Device you want to print the Install messages.
You can queue the install by enter a 'Q' at the device prompt.
Enter a '^' to abort the install.

DEVICE: HOME// ;;99999  HOME  (CRT)
-----

Install Started for MOCHA ENH 2 COMBINED BUILD 1.0 :
      Aug 11, 2015@15:07:15

Build Distribution Date: Aug 11, 2015

Installing Routines:
      Aug 11, 2015@15:07:15

Install Started for GMRA*4.0*46 :
      Aug 11, 2015@15:07:15

Build Distribution Date: Aug 11, 2015

Installing Routines:
      Aug 11, 2015@15:07:15

Running Post-Install Routine: POST^GMRAY46

Updating Routine file...

Updating KIDS files...

GMRA*4.0*46 Installed.
      Aug 11, 2015@15:07:15

Not a production UCI

NO Install Message sent

Install Started for OR*3.0*269 :

```

```

Aug 11, 2015@15:07:15

Build Distribution Date: Aug 11, 2015

Installing Routines:
    Aug 11, 2015@15:07:15

Running Pre-Install Routine: PRE^ORY269

Deleting existing ORDER CHECK INSTANCES data dictionary while preserving data...
DONE

Installing Data Dictionaries:
    Aug 11, 2015@15:07:15

Running Post-Install Routine: POST^ORY269 .

Order Check Expert System Rule Transporter
Created: JUL 26,2013 at 06:56 at CPRS30.FO-SLC.MED.VA.GOV
Current Date: AUG 11,2015 at 15:07 at CHEYL37.FO-BAYPINES.MED.VA.GOV

Loading Data . . . . .
  Installing '863.8 OCX MDD PARAMETER' records... . . . . .
  . . . . .
  Installing '864.1 OCX MDD DATATYPE' records... . .
  Installing '863.4 OCX MDD ATTRIBUTE' records... . . . .
  Installing '863.2 OCX MDD SUBJECT' records... .
  Installing '863.3 OCX MDD LINK' records... . . . .
  Installing '860.9 ORDER CHECK NATIONAL TERM' records... . . . .
  . . . . .
  Installing '860.8 ORDER CHECK COMPILER FUNCTIONS' records... . . . .
  . . . . .
  Installing '860.6 ORDER CHECK DATA CONTEXT' records... . . .
  Installing '860.5 ORDER CHECK DATA SOURCE' records... . . .
  Installing '860.4 ORDER CHECK DATA FIELD' records... .

    No data filing errors.

Transport Finished...

---Creating Order Check Routines-----

Build list of Active Rules, Elements and Datafields...

    96 DATA FIELDS

    81 ELEMENTS

    39 RULES

Compile DataField Navigation code...

    101 DataField Navigation Code Arrays

Compile Element Evaluation code...

    75 Event Evaluation Code Arrays

Compile Element MetaCode...

    81 Element Metacode Arrays

```

```

Get Compiler Function Code...

    51 Compiler Include Functions
Compile Rule Element Relation code...

    56 Rule Element Relation Code Arrays
Construct Decision Tree...

    716 Sub-Routines
Optimize Sub-Routines...

    287 Sub-Routines

    60% Optimization
Assemble Routines...

    39 OCXOZ* Routines
Converting existing data in the ORDER CHECK INSTANCES file...
DONE

Updating Routine file...

Updating KIDS files...

OR*3.0*269 Installed.
    Aug 11, 2015@15:07:26

Not a production UCI

NO Install Message sent

Install Started for PSJ*5.0*281 :
    Aug 11, 2015@15:07:26
Build Distribution Date: Aug 11, 2015

Installing Routines:
    Aug 11, 2015@15:07:26

Installing PACKAGE COMPONENTS:

Installing PROTOCOL
    Aug 11, 2015@15:07:26

Updating Routine file...

Updating KIDS files...

PSJ*5.0*281 Installed.
    Aug 11, 2015@15:07:26

Not a production UCI

NO Install Message sent

Install Started for PSO*7.0*411 :
    Aug 11, 2015@15:07:26

```

```
Build Distribution Date: Aug 11, 2015

Installing Routines:
    Aug 11, 2015@15:07:26

Running Post-Install Routine: ADDCR^PSOCROC

Updating Routine file...

Updating KIDS files...

PSO*7.0*411 Installed.
    Aug 11, 2015@15:07:26

Not a production UCI

NO Install Message sent

Updating Routine file...

Updating KIDS files...

MOCHA ENH 2 COMBINED BUILD 1.0 Installed.
    Aug 11, 2015@15:07:26

No link to PACKAGE file

Install Completed
```

## 7.3 Routine Deletions

After the patch is successfully installed, KIDS will automatically delete the main post-install routine ORY269. Sites may safely delete the following remaining post-install routines using the Delete Routines [XTRDEL] option:

```
ORY2690
ORY26901
ORY26902
ORY26903
ORY26904
ORY26905
ORY2691
ORY2692
ORY2693
ORY2694
ORY269ES
```

## 8. Rollback Procedure

This last section of the document is only needed if there are major problems resulting from the installation of these patches, and there is concurrence from the Enterprise Product Support and MOCHA Development Teams that a rollback must occur.

If you are a CMOP (Consolidated Mail Outpatient Pharmacy), and the only build you installed was the PSS\*1.0\*175 build, all you will need to do is restore your routine backup made of that install. Then you can delete the routines PSSCKOS and PSSDIUTX, which were new with PSS\*1.0\*175, and you are done.

If a rollback occurs, entries created in the ORDER CHECK INSTANCES file (#100.05) while ME2 (MOCHA Enhancements 2) was installed are retained. However, data in fields that are new with ME2 is lost and unrecoverable.

The CLINICAL REMINDER entry that was added to the APSP INTERVENTION TYPE file (#9009032.3) in the PSO\*7.0\*411 Post-Init routine will remain in the file.

### 8.1 Rollback Requirements

This section outlines the requirements for rollback.

#### 8.1.1 Software Required for Rollback

Make sure you have the following software available for the rollback:

1. Routine backup from PSS\*1.0\*175 Install
2. Routine backup from MOCHA\_ENH\_2\_GMRA\_OR\_P SJ\_PSO install
3. MOCHA\_ENH\_2\_ROLLBACK.KID Host File
4. (Local) Your local backups of the ORDER CHECK INSTANCES (#100.05) File Data Dictionary and the PSJ DISPLAY DRUG ALLERGIES Protocol from Section 5, if a backup was made for local modifications.

\*\*\*Note – The MOCHA\_ENH\_2\_ROLLBACK.KID build is available from Enterprise Product Support upon request. Enterprise Product Support can also supply the routine backups if needed.

#### 8.1.2 Overview of Rollback Build

The MOCHA\_ENH\_2\_ROLLBACK combined build is made up of the following single builds:

1. ORZ\*3.0\*269
  - Requires OR\*3.0\*352, OR\*3.0\*306, OR\*3.0\*311, OR\*3.0\*346, OR\*3.0\*286, GMRA\*4.0\*46, OR\*3.0\*395
  - ORDER CHECK INSTANCES (#100.05) File Data Dictionary restores the data dictionary for the ORDER CHECK INSTANCES file (#100.05) back to its pre-OR\*3\*269 state.
  - Converts existing data in the ORDER CHECK INSTANCES file (#100.05) into a pre-OR\*3\*269 format. Data in new fields is deleted.
  - Restores the ALLERGY - CONTRAST MEDIA REACTION order check back to its pre-OR\*3\*269 state.
  - ORY269 conversion routines that will be deleted after patch installation
  - Instruction provided in section 8.6 of this guide on other routines to delete

2. PSJZ\*5.0\*281

- Requires PSJ\*5.0\*281
- PSJ DISPLAY DRUG ALLERGIES Protocol to restore the Protocol without the EXIT ACTION Field data.
- Deletes the following routines that were new with the Mocha Enhancements 2 install:
  - o PSJCROC
  - o PSJNEWOA
  - o PSJNEWOC

3. PSOZ\*7.0\*411

- Requires PSO\*7.0\*411
- Deletes the following routines that were new with the Mocha Enhancements 2 install:
  - o PSOCROC
  - o PSODGAL3
  - o PSONEWOA
  - o PSONEWOC
  - o PSOOCKV1

4. PSSZ\*1.0\*175

- Requires PSS\*1.0\*175
- Deletes the following routines that were new with the PSS\*1.0\*175 install:
  - o PSSCKOS
  - o PSSDIUTX

## 8.2 Rollback Considerations/Restrictions

The rollback should occur when Pharmacy applications are not in use, no other pharmacy patches are being installed and when tasked jobs from Clinical Applications are not running.

Rollback should also occur when CPRS usage is at a minimum, particularly Pharmacy activities. Rollback should occur when allergy data is not being entered or edited in CPRS.

The length of time for the rollback is dependent upon the number of entries in the ORDER CHECK INSTANCES (#100.05) file that will be converted and could range from five to thirty minutes.

**Do not queue any of the rollback process.**

Failure to heed this warning may result in 'source routine edited' errors during a database update. Edits may be lost and records may be left in an inconsistent state.

## 8.3 Rollback Steps


1. Restore the routine backup from the patch PSS\*1.0\*175 install.
2. Restore the routine backup from the MOCHA\_ENH\_2\_GMRA\_OR\_PSJ\_PSO install.

\*\*\*Note - if for some reason the MOCHA\_ENH\_2\_GMRA\_OR\_PSI\_PSO build was installed multiple times, make sure you restore the back up from the first install.

3. Use **Load a Distribution**. You may need to prepend a directory name. When prompted for "Enter a Host File:", respond with MOCHA\_ENH\_2\_ROLLBACK.KID.

Example: USER\$:[ABC]MOCHA\_ENH\_2\_ROLLBACK.KID

4. From this menu, you may then select to use the following options:

<b>Note:</b>	The following are OPTIONAL - (When prompted for the INSTALL NAME, enter MOCHA_ENH_2_ROLLBACK 1.0).
	

- **Verify Checksums in Transport Global** - This option will allow you to ensure the integrity of the routines that are in the transport global.
  - **Print Transport Global** - This option will allow you to view the components of the KIDS build.
  - **Compare Transport Global to Current System** - This option will allow you to view all changes that will be made when this patch is installed. It compares all components of this patch (routines, DD's, templates, etc.).
  - **Backup a Transport Global** - This option will create a backup message of any routines exported with this patch. It will not backup any other changes such as DDs or templates.
5. Select the Installation option Install Package(s). This is the step to start the installation of this KIDS patch:
    - a. Choose the Install Package(s) option to start the patch install and enter MOCHA\_ENH\_2\_ROLLBACK 1.0 at the INSTALL NAME prompt.
    - b. When prompted 'Want KIDS to Rebuild Menu Trees Upon Completion of Install? NO//', answer NO.
    - c. When prompted 'Want KIDS to INHIBIT LOGONs during the install? NO//' answer NO.
    - d. When prompted 'Want to DISABLE Scheduled Options, Menu Options, and Protocols? NO//' answer YES.
    - e. When prompted 'Enter options you wish to mark as 'Out Of Order':' answer GMRA PATIENT A/AR EDIT.
    - f. When prompted 'Enter protocols you wish to mark as 'Out Of Order':' press the Enter key.
    - g. When prompted 'Delay Install (Minutes): (0-60): 0//' press the Enter key.
  6. Local Modifications (If Applicable):

Install your local backups of the ORDER CHECK INSTANCES (#100.05) File Data Dictionary and the PSJ DISPLAY DRUG ALLERGIES Protocol from Section 5, if a backup was made.

If either of the 2 routine backups (PSS\*1.0\*175 or MOCHA\_ENH\_2\_GMRA\_OR\_PSI\_PSO) had to be retrieved from Enterprise Product Support, any local modifications can now be re-entered.

#### 7. Patch Installation Tracking:

The site will need to work with the Enterprise Product Support and MOCHA Development Teams to make the necessary changes to the PACKAGE (#9.4) File to accurately reflect the current state of the account in regards to the patches installed.

The MOCHA Development Team will provide instructions for this step to Enterprise Product Support in the 'ME2 Rollback - Update Patch Application Tracking' document.

#### 8. Routine Deletions:

After the patch is successfully installed, KIDS will automatically delete the main post-install routine ORY269. Sites may safely delete the following remaining post-install routines using the Delete Routines [XTRDEL] option:

ORY2690  
ORY26901  
ORY26902  
ORY26903  
ORY26904  
ORY26905  
ORY26906  
ORY26907  
ORY2691  
ORY2692  
ORY2693  
ORY2694  
ORY269ES

Note: If the Kernel parameter XPD NO\_EPP\_DELETE is set to 1 in the account, the routines ORY269, ORZY269, ORZY269E and GMRAY46 will *not* be deleted by KIDS during the installation. For purposes of rollback however; they can be deleted using this option.



## Appendix: Acronyms

Term	Definition
CPRS	Computerized Patient Record System
DD	Data Dictionaries
FDB	First Databank
FTP	File Transfer Protocol
HWSC	Heath <u>e</u> Vet Web Services Client
IRMS	Information Resources Management Service
KIDS	Kernel Installation and Distribution System
MOCHA	Medication Order Check Healthcare Application
MUMPS	Massachusetts General Hospital Utility Multi-Programming System
NDF	National Drug File
PDM	Pharmacy Data Management
PECS	Pharmacy Enterprise Customization System
PEPS	Pharmacy Enterprise Product System
PMI	Patient Medication Information
VA	Department of Veterans Affairs
VDL	VA Software Document Library
VISTA	Veterans Health Information Systems and Technology Architecture