



Update_2_0_91

Clinical Reminders

VA-VVC BLOOD PRESSURE TEMPLATE

Install Guide

October 2019

Product Development
Office of Information Technology
Department of Veterans Affairs

Table of Contents

INTRODUCTION	1
INSTALL DETAILS	2
INSTALL EXAMPLE.....	2
POST INSTALLATION.....	6

Introduction

The VA Video Connect Blood Pressure (VVC BP) template will be used primarily by PACT nursing staff to complete virtual BP check visits with Veterans. It may also be used by providers (PACT, nephrology, cardiology, etc or by Clinical Pharmacy Specialists) occasionally.

The process for completing a virtual BP check may involve issuing a new BP cuff to the Veteran via mail, instructing the Veteran on proper usage, measuring the BP more than once (once in each arm if the initial visit), medication review, symptom review, and education regarding behavior modifications for optimal BP management.

This new process for BP monitoring now allows clinicians to record a remotely measured BP into the vital signs package. This was reviewed and approved by the Performance Accountability Workgroup in May 2018.

Veterans will decide with their clinician whether they are appropriate for and prefer a virtual visit vs an in-person visit.

The template will be connected to a note title (VA VIDEO CONNECT BLOOD PRESSURE VISIT). **At this time no health factors are being collected.**

UPDATE_2_0_91 contains 1 Reminder Exchange entry:

UPDATE_2_0_91 VA-VCC BLOOD PRESSURE

The exchange file contains the following components:

TIU TEMPLATE FIELD

BLANK LINE 1
WORD
BLANK SP
OTHER TEXT

HEALTH FACTORS

VA-REMINDER UPDATES
VA-UPDATE_2_0_91

REMINDER SPONSOR

Office of Connected Care

REMINDER TERM

VA-REMINDER UPDATE_2_0_91

REMINDER DIALOG

VA-TH VVC HYPERTENSION

Install Details

This update is being distributed as a web host file. The address for the host file is:
http://vista.med.va.gov/reminders/UPDATE_2_0_91.PRD

The file will be installed using Reminder Exchange, programmer access is not required.

Installation:

=====


This update can be loaded with users on the system. Installation will take less than 15 minutes.

Install Example

If your site is unable to use LWH functionality, you will use vista mail to install the exchange file.

To Load the Web Host File. Navigate to Reminder exchange in Vista

```
+          + Next Screen  - Prev Screen  ?? More Actions  >>>
CFE  Create Exchange File Entry          LHF  Load Host File
CHF  Create Host File                    LMM  Load MailMan Message
CMM  Create MailMan Message              LWH  Load Web Host File
DFE  Delete Exchange File Entry          RI   Reminder Definition Inquiry
IFE  Install Exchange File Entry          RP   Repack
IH   Installation History
Select Action: Next Screen// LWH  Load Web Host File
Input the URL for the .prd file: http://vista.med.va.gov/reminders/UPDATE_2_0_91.PRD
```



At the **Select Action:** prompt, enter **LWH** for Load Web Host File

At the **Input the url for the .prd file:** prompt, type the following web address:

http://vista.med.va.gov/reminders/UPDATE_2_0_91.PRD

```
http://vista.med.va.gov/reminders/UPDATE_2_0_91.PRD successfully loaded.
```

You should see a message at the top of your screen that the file successfully loaded.
Search and locate an entry titled **UPDATE_2_0_91 VA-VCC BLOOD PRESSURE** in reminder exchange

+Item	Entry	Source	Date Packed
184	UPDATE_2_0_91 VA-VCC BLOOD PRESSURE	GRIFFITH@SALT LAKE CI	09/17/2019@14:03
185	UPDATE_2_0_92 VA-BFA TEMPLATE UPDATE	GRIFFITH@SALT LAKE CI	08/06/2019@07:14
186	UPDATE_2_0_93 VA-VIONE POLYPHARMACY REVIEW	PLOTT@SALT LAKE CI	09/05/2019@12:22
187	UPDATE_2_0_94 VA-CRC UPDATE	GRIFFITH@SALT LAKE CI	08/28/2019@11:50
188	UPDATE_2_0_95 VA-COMMUNITY CARE EMERGENCY CARE UPDATE	GRIFFITH@SALT LAKE CI	08/28/2019@13:09
189	UPDATE_2_0_96 VA-CAREGIVER	GRIFFITH@SALT LAKE CI	09/06/2019@13:09
+ + Next Screen - Prev Screen ?? More Actions >>>			
CFE	Create Exchange File Entry	LHF	Load Host File
CHF	Create Host File	LMM	Load MailMan Message
CMM	Create MailMan Message	LR	List Reminder Definitions
DFE	Delete Exchange File Entry	LWH	Load Web Host File
IFE	Install Exchange File Entry	RI	Reminder Definition Inquiry
IH	Installation History	RP	Repack
Select Action: Next Screen// IFE Install Exchange File Entry			
Enter a list or range of numbers (1-297): 184			

At the **Select Action** prompt, enter **IFE** for Install Exchange File Entry
Enter the number that corresponds with your entry titled **UPDATE_2_0_91 VA-VCC BLOOD PRESSURE** (in this example it is entry 184, it will vary by site)

Component	Category	Exists
Source: GRIFFITH,ELIZABETH A at SALT LAKE CITY		
Date Packed: 09/17/2019@14:03:25		
Package Version: 2.0P35		
Description:		
The following Clinical Reminder items were selected for packing:		
REMINDER DIALOG		
VA-TH VVC HYPERTENSION		
REMINDER TERM		
VA-REMINDER UPDATE_2_0_91		
Non-exchangeable order dialog(s):		

Name: GMRAOR ALLERGY ENTER/EDIT		
Type: Action		
Display Text: Allergy/Adverse Reaction		
+ Enter ?? for more actions >>>		
IA	Install all Components	IS Install Selected Component
Select Action: Next Screen// IA		

At the **Select Action** prompt, type **IA** for Install all Components and hit enter.
Select Action: Next Screen// **IA Install all Components**

During the installation, you will see several prompts: *Note that you will **Install** all new components*

You will be prompted to install the dialog component:

Packed reminder dialog: VA-TH VVC HYPERTENSION [NATIONAL DIALOG]				
Item	Seq.	Dialog Findings	Type	Exists
1		VA-TH VVC HYPERTENSION	dialog	
2	5	VA-VCC NOTE VERSION BP Finding: *NONE*	element	
3	10	VA-GP TH VVC HYPERTENSION DISCLOSURE Finding: *NONE*	group	
4	10.5	VA-GP TH VVC HYPERTENSION LOCATION Finding: *NONE*	group	
5	10.5.5	VAL-TH VVC HYPERTENSION LOCATION HOME Finding: *NONE*	element	
6	10.5.10	VA-TH VVC HYPERTENSION LOCATION OTHER Finding: *NONE*	element	
7		PXRM COMMENT	prompt	X
8	10.10	VA-GP TH VVC HYPERTENSION EMERGENCY	group	
+ + Next Screen - Prev Screen ?? More Actions				
DD	Dialog Details	DT	Dialog Text	IS Install Selected
DF	Dialog Findings	DU	Dialog Usage	QU Quit
DS	Dialog Summary	IA	Install All	
Select Action: Next Screen// IA				

At the **Select Action** prompt, type **IA** to install the dialog – **VA-TH VVC HYPERTENSION**

Select Action: Next Screen// **IA Install All**

Install reminder dialog and all components with no further changes: **Y// ES**

Packed reminder dialog: VA-TH VVC HYPERTENSION [NATIONAL DIALOG]				
VA-TH VVC HYPERTENSION (reminder dialog) installed from exchange file.				
Item	Seq.	Dialog Findings	Type	Exists
1		VA-TH VVC HYPERTENSION	dialog	X
2	5	VA-VCC NOTE VERSION BP Finding: *NONE*	element	X
3	10	VA-GP TH VVC HYPERTENSION DISCLOSURE Finding: *NONE*	group	X
4	10.5	VA-GP TH VVC HYPERTENSION LOCATION Finding: *NONE*	group	X
5	10.5.5	VAL-TH VVC HYPERTENSION LOCATION HOME Finding: *NONE*	element	X
6	10.5.10	VA-TH VVC HYPERTENSION LOCATION OTHER Finding: *NONE*	element	X
7		PXRM COMMENT	prompt	X
8	10.10	VA-GP TH VVC HYPERTENSION EMERGENCY	group	X
+ + Next Screen - Prev Screen ?? More Actions				
DD	Dialog Details	DT	Dialog Text	IS Install Selected
DF	Dialog Findings	DU	Dialog Usage	QU Quit
DS	Dialog Summary	IA	Install All	
Select Action: Next Screen// Q				

When the dialog has completed installation, you will then be returned to this screen. At the

Select Action prompt, type **Q**.

Component	Category	Exists
Source: GRIFFITH,ELIZABETH A at SALT LAKE CITY		
Date Packed: 09/17/2019@14:03:25		
Package Version: 2.0P35		
Description:		
The following Clinical Reminder items were selected for packing:		
REMINDER DIALOG		
VA-TH VVC HYPERTENSION		
REMINDER TERM		
VA-REMINDER UPDATE_2_0_91		
Non-exchangeable order dialog(s):		

Name: GMRAOR ALLERGY ENTER/EDIT		
Type: Action		
Display Text: Allergy/Adverse Reaction		
+ Enter ?? for more actions >>>		
IA	Install all Components	IS Install Selected Component
Select Action: Next Screen// Q		

You will be returned to this screen. Type Q to quit.

Install Completed

Post Installation

1. Make the Dialog templates available to be attached to a Personal or Shared template or to be able to attach the dialog to a progress note title.

From the Reminders Manager Menu:

Select Reminder Managers Menu Option: **CP CPRS Reminder Configuration**

CA Add/Edit Reminder Categories
CL CPRS Lookup Categories
CS CPRS Cover Sheet Reminder List
MH Mental Health Dialogs Active
PN Progress Note Headers
RA Reminder GUI Resolution Active
TIU TIU Template Reminder Dialog Parameter
DL Default Outside Location
PT Position Reminder Text at Cursor

Select CPRS Reminder Configuration Option: **TIU TIU Template Reminder Dialog Parameter**

Reminder Dialogs allowed as Templates may be set for the following:

- 1 User USR [choose from NEW PERSON]
- 3 Service SRV [choose from SERVICE/SECTION]
- 4 Division DIV [choose from INSTITUTION]
- 5 System SYS [NATREM.FO-SLC.MED.VA.GOV]

Enter selection: **5** *Choose the appropriate number for your site. Your site may do this by System or other levels. The example below uses SYSTEM level*

Setting Reminder Dialogs allowed as Templates for System:

Select Display Sequence: **?**

148 VA-VETERANS CHOICE NOTE
149 VA-TBI SCREENING
151 VACO TRANSPLANT REFERRAL (D)

*When you type a question mark above, you will see the list of #'s (with dialogs) that are already taken. Choose a number **NOT** on this list. For this example, looking above I see the number 150 is not present, so I will use 150.*

Select Display Sequence: **150**

Are you adding 150 as a new Display Sequence? Yes// **Y YES**

Display Sequence: 150// **<Enter>**

Clinical Reminder Dialog: **VA-TH VVC HYPERTENSION Enter>**

OK? Yes// **<Enter>** (Yes)

2. Setup of Note title

The new note title that needs to be created is:

VIDEO BLOOD PRESSURE VISIT

This note title is to be mapped to the following VHA Enterprise Standard Title:
TELEHEALTH NOTE

Comprehensive information on Creating Document Definitions can be found beginning on page 47 of the [TIU/ASU Implementation Guide](#)

3. Associate the reminder dialog with the note title in CPRS


To do this, go to Template Editor, then Document Titles, select the appropriate folder (if applicable), and choose New Template (top right-hand corner).

Type in the following information:

Name: VIDEO BLOOD PRESSURE VISIT

Template Type: Reminder Dialog

Reminder Dialog: VA-TH VVC HYPERTENSION



Next, type in the Associated Title, which in this case is VIDEO BLOOD PRESSURE VISIT and hit Apply.

4. Add your local DATA OBJECT FOR PATIENT ADDRESS to the following element:
VAL-TH VVC HYPERTENSION LOCATION HOME
Using the menu Reminder Dialog Management, Reminder Dialogs, change your view to Element.
Search for and edit the following element: VAL-TH VVC HYPERTENSION LOCATION HOME.
5. Add your local DATA OBJECT FOR PATIENT PHONE NUMBER to the following element: VAL-TH VVC HYPERTENSION PHONE NUMBER OBJECT
Using the menu Reminder Dialog Management, Reminder Dialogs, change your view to Element.
Search for and edit the following element: VAL-TH VVC HYPERTENSION PHONE NUMBER OBJECT
6. Add your local DATA OBJECTS FOR MEDICATIONS AND ALLERGIES to the following element: VAL-TH VVC HYPERTENSION DISCLOSURE MED AND ALLERGY OBJECT
Using the menu Reminder Dialog Management, Reminder Dialogs, change your view to Element.
Search for and edit the following element: VAL-TH VVC HYPERTENSION DISCLOSURE MED AND ALLERGY OBJECT
7. Add a local DATA OBJECT FOR BMP WITHIN THE LAST 2 YEARS to the following element: VAL-TH VVC HYPERTENSION DISCLOSURE LAB OBJECT
Using the menu Reminder Dialog Management, Reminder Dialogs, change your view to Element.
Search for and edit the following element: VAL-TH VVC HYPERTENSION DISCLOSURE LAB OBJECT

NOTE: If you require further technical assistance, if there are any questions or problems during the installation, please notify your local IT support to log a national help desk ticket Desk and have them submit a national ticket to:

Category: Enterprise Applications

Enterprise Application: VistA - Clinical Reminders

Assignment Group: NTL SUP Clin 2